



## DYNAMIC LEG EXTENSION DEVICE

This device was designed for users who display strong extensor spasticity leading to ongoing breakage of wheelchair frames or standard footrest hangers and mounts. This device absorbs and diffuses tone by providing progressive resistance during the downward movement of the lower extremity and then returns the lower extremity to a preferred position.

- Provides progressive resistance from 5 lbs – 60 lbs (Pediatric) or 17 lbs – 120 lbs (Adult)
- Allows for up to 2.75" of Dynamic Leg Extension
- Mounts to 1" OD footrest hanger tubes
- Direct receiver mount for Miller's Leg Extension Tube models #33150, #33151, #33153 & #33157
- May fit other manufacturer's leg extension tubes if they are 3/4" OD

### Miller's Dynamic Leg Extension Devices

HCPC: K0108

- **33600** (Dyn Leg-Ext Dev, Ea)
- **33600H** (Dyn Leg-Ext Dev, High Mnt, Ea)
- **33600P** (Dyn Leg-Ext Ped Dev, Ea)
- **33600PH** (Dyn Leg-Ext Ped Dev, High Mnt, Ea)

Please contact our sales staff to assist you in placing an order for any of our products or to answer any questions that you may have. If you can't find what you need in our product catalog or on our website, we can work with you to create a custom solution.



*We specialize in creating solutions to expand your capabilities!*

2023 Romig Rd. Akron, OH 44320  
Toll Free (800)837-4544  
Fax (866)572-2603  
**MillersAdaptive.com**